## smart HCY Test

Testkit for quantitative determination of L-Homocysteine in serum on smart Photometer.

**EUROLyser**  
DIAGNOSTICA

Eurolyser Diagnostica GmbH Tel. +43 (0) 662 43 21 00  
Bayernstraße 11a Fax. +43 (0) 662 43 21 00 50  
5020 Salzburg/Austria www.eurolyser.com

Order Information	Indication	Kit size
Order No.: ST0120	smart HCY-Test Kit	16 tests/kit
Order No.: ST1200	smart HCY-Control Kit	2 x 1 ml / high and low
Order No.: ST1210	smart HCY-Calibrator Kit	1 x 1 ml

### English:

#### Test kit preparation:

Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put back test kit package into refrigerator.

### Summary:

Homocysteine is a thiol containing Amino Acid produced by the intracellular demethylation of methionine. Total Homocysteine tHCY represents the sum of all forms of HCY including forms of oxidized, protein bound and free. Elevated level of tHCY has emerged as an important risk factor in the assessment of cardiovascular disease. Excess HCY in the blood stream may cause injuries to arterial vessels due to its irritant nature and result in inflammation and plaque formation, which may eventually cause blockage of blood flow to the heart. Elevated tHCY levels are caused by four major factors including:

- genetic deficiencies
  - nutritional deficiency in B vitamins such as B6 B12 and folate
  - renal failure for effective amino acid clearance
  - drug interactions such as nitricide, methotrexate and phenytoin
- Elevated levels of tHCY are also linked with Alzheimer's disease and Osteoporosis

### Assay principle:

Kinetic determination based on enzymatic Cycling Methode:  
Serum: 3-50µmol/L (small lot specific variations possible)

### Testkit:

ERS Cuvette filled with 1000µl Co-Substrate Sam 0,1mM  
ERS Cap filled with 250µl Enzymreagent with SAH Ado glutamate dehydrogenase

### Stability and Storage:

Store at 2-8°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels

### Warning:

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds.

### Sample Material:

Fresh serum is recommended for the HCY assay. It is important to centrifuge samples immediately after collection to separate the serum from the blood cells. Hemolysed or turbid specimens or severely lipemic specimens are not recommended for this HCY assay.

### Reference Ranges:

adults < 12µmol/L ; Each laboratory is recommended to establish a range of normal values for the population in their region.

### Control Material:

A set of normal and abnormal HCY controls is available: HCY Control Kit Order No.: ST1200

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