## smart CRP Test

In Vitro Diagnostic test kit for determination of c-reactive protein in serum or plasma and whole blood with smart photometer.

### Order Information

Order No.: ST0100  
Order No.: ST1000  
Order No.: ST1010

### Indication

smart CRP-Test Kit  
smart CRP-Control Kit  
smart CRP-Calibrator Kit

### Kit size

32 tests/pack  
2 x 1 ml (low high range)  
1 x 1 ml

### English:

#### Test kit preparation:

Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

### Summary:

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. It increases already after 6 hours in acute inflammatory processes associated with bacterial infections, post operative conditions or tissue damage. Many studies have demonstrated that in apparently healthy subjects there is a direct correlation between CRP concentrations and the risk for developing coronary heart disease (CHD).

### Methode-Measurement Range:

Immunturbidimetric Assay:

Serum/Plasma: 0,05-12 mg/dl (0,5-120mg/l) for SMART546 and 0,10-12 mg/dl (1,0-120mg/l) for SMART700  
whole blood: 0,20-20 mg/dl (2,0-200mg/l) for SMART546 and 0,20-24 mg/dl (2,0-240mg/l) for SMART700

### Principle:

Kinetic determination of the concentration of CRP by photometric measurement at 546nm or 700nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

### Test Kit:

ERS Cuvette pre-filled with 1000µl buffer reagent Glycine Buffer 170mM  
ERS Cap pre-filled with 250µl latex reagent with rabbit anti human CRP antibodies (0,2%)

### Storage Instructions and Stability of the Test Kit:

The reagents are stable until the indicated month of expiry , if stored at 2-8°C.

### Warnings and Precautions:

The reagents contain sodium azide (0,95g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.. Take the necessary precautions for the use of laboratory reagents.

### Waste Management:

please refer to local legal requirements.

### Specimen Sample Material:

Freshly drawn serum is preferred.  
Heparin plasma, EDTA plasma, edta whole blood and capillary blood from finger tip may also be used if limited test performance is accepted.

### Reference Range:

adults < 5mg/l (<0,5 mg/dl) ; values > 3mg/l (> 0,3 mg/dl) can indicate a risk for coronary heart disease CHD.  
It is recommended that each laboratory establish its own expected range.

### Quality Control:

for internal quality control the smart CRP Control Calibrator Kit order no.:ST1000 should be used

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### 2.) Attention: Select whole blood or serum via the analyser menu

- For detailed information please refer to the user manual of the smart laboratory photometer
- 2.1) Sample type menu
  - 2.2) Hematokrit
  - 2.3) Starting measurement

### 3.) Sample preparation for whole blood (fingertip-or EDTA primary tube)

- 3.1a) Aspirate 5µl from fingertip OR ...
- 3.1b) Aspirate 5µl from EDTA capillary blood collection system OR ...
- 3.1c) Aspirate 5µl from EDTA primary tube
- 3.2) Dispense 5µl sample into liquid in ERS cuvette
- 3.3) Apply ERS cap

### ATTENTION!

Do not wipe the pipette tip after aspirating the sample!  
Incubate tests 10 minutes at roomtemperature before use!

### 4.) Sample preparation for Serum or Plasma (primary tube serum , Li-Hep or EDTA Plasma)

- 4.1a) Aspirate 5µl serum from primary tube OR ...
- 4.1b) Aspirate 5µl plasma from centrifuged capillary blood collection system
- 4.2) Dispense 5µl sample into liquid in ERS cuvette
- 4.3) Apply ERS cap

### 5.) Sample processing with smart photometer

Please follow the instructions in the users manual